

# iOS Update



1. Plug your device into power and connect to the internet with WiFi.
2. Go to Settings > General, then tap Software Update.



3. Tap Download and Install
4. To update now, tap Install.

Or, you can tap Later and choose: Install Tonight or Remind Me Later.

If you tap Install Tonight, just plug your device into power before you go to sleep. Your device will update automatically overnight.

## Important

Some clinic VPN or proxy connections might prevent your device from contacting the update servers. In this case - please use a WiFi connection with-out such restrictions.

Need more help? Contact us on 1300 535 645 or [contactus@ortesi.com.au](mailto:contactus@ortesi.com.au)

**ortesi**  
3D printed braces & orthotics

Ortesi Pty Ltd  
ABN 14 153 078 936  
P.O Box 4555  
Dandenong South VIC 3164  
P 1300 535 645  
F 1300 515 685  
[contactus@ortesi.com.au](mailto:contactus@ortesi.com.au)