

# LABORATORY GUIDE MODIFICATIONS



MODIFICATIONS	MILD	MODERATE	AGGRESIVE
MEDIAL SKIVE	2mm	3mm	4mm
LATERAL SKIVE	2mm	3mm	4mm
PLANTAR FASCIAL ACCOM.	1mm	2mm	4mm
CUBIOD NOTCH	2mm	3mm	4mm
PLANTER 5 <sup>TH</sup> RAY GRIND (AKA LATERAL CAST GRIND)	4mm	5mm	6mm
LABORATORY STANDARDS	LOW	STANDARD	HIGH
HEEL CUP HEIGHT (PP)	6mm	12mm	18mm
HEEL CUP HEIGHT (EVA)	4mm	8mm	14mm