


nora Lunasoft SL | ortesi

COVER RANGE

3D printed braces & orthotics

Genuine Nora Lunasoft SL is an extremely lightweight, smooth, resilient and dimensionally stable covering material. Thanks to its closed cell structure it is washable, extremely hygienic and skin friendly. Unlike standard EVA cover materials, genuine Nora Lunasoft SL will not bottom out over time.

The range has a standard thickness of 2.0mm, shore hardness of A40, and density of 0,20g/cm³.



BLACK



BEIGE



GREY



BLUE



RED



GREEN



DARK BLUE




ORANGE



BLACKBERRY




LIME



DARK BROWN



RASPBERRY



CAMOUFLAGE



QAUTRO



CORK



BROWN & LIME




ORANGE & BLUE




BLACK & WHITE




BLACK & LILAC




BLACK & PINK




BLACK & YELLOW



BLACK & RED



SAFARI



TRI-COLOUR