


PREMIUM COVER RANGE

ortesi
3D printed braces & orthotics

NATURAL LEATHER – SOFT

Our NATURAL LEATHER – SOFT is sourced from South Korea. The hides are carefully procured from Heifer cows after which the hides are cold water treated after dyeing to enhance the grain in increase softness.




BLACK



BROWN

NATURAL LEATHER – FIRM

Our NATURAL LEATHER – FIRM is sourced from Pakistan. The heavy nature of the material results in higher abrasion resistance thanks to its heavy pigmentation.



BLACK




BEIGE

SYNTHETIC LEATHER

Our SYNTHETIC LEATHER is an **ECO FRIENDLY** man made leather.


The material is fully bio-degradable and pre-treated with GREENGAURD stain resistant barrier. Made from 85% H.D.P.U and 15% Polyester.



BLACK

P.S. VLIES


P.S. VLIES is a hard-wearing, non-woven cotton-polyester polyurethane blend. The material absorbs moisture during wearing and dries effectively afterwards without any negative effects.



BLACK

SPENCO 1.5mm & 3.0mm

SPENCO is an expanded black neoprene with a Nylon stretch cover adhered to one side. The material has excellent cushioning and impact properties and is super-soft to the touch. The material is shock absorbing and extremely comfortable making it an ideal choice for active patients.



BLACK



BLUE



GREEN

NOTE – GREEN IS AVAILABLE IN 3.0mm ONLY